

## PLANTAR WARTS

For the treatment of warts, K-Laser's objective is to propose a non-invasive therapy solution that consists in a single, or double treatment cycles, depending on wart's characteristics.



In the images a female patient affected by a wart on the feet, previously treated with salicylic acid.

The warts were coloured with a dark pigment before the laser applications.

The protocol selected was 'Wart adult' in the Podo-dermatology area with the following modality: no-contact surgical handpiece directly on the wart with a peak power of 20W, T-on 18 ms, T-off 54 ms, pulsed modality, average power of 4.9 W and total delivered energy of 360 mJ.